

BARIATRIC TRANSFER MATRIX

			TARGET SURFACE measure			
			SUPINE postures	SITTING postures	STANDING postures	
Functional Independence Measure	STARTING SURFACE postures	Supervision & Verbal Clues – Min, Mod, Max & Total Assist	SUPINE	Slide Glide or Sling Devices, Transfer Sheets, Transfer Mattress, Roller Boards, Supine Plastic Transfer Boards, Trapeze	Sling Devices: Bed Side Lifts, A-Frame Lifts, Trapeze, Other Sling Devices	Stand Frame or Slide Glide to Tilt Table
			SITTING	Sling Devices	Sling Devices	Standing Frame
			FLOOR	Sling Devices & Pneumatic Risers	Sling Devices	Not Appropriate
			SUPINE	Guarded Supine Self-Scooting	Guarded Sitting Scooting or Stand Pivot with Gait Belt	Guarded Stand Pivot Transfer Safety Gait Belt and Baseline Gait Device
			SITTING	Guarded Sitting Self-Scooting	Guarded Sitting Scooting or Stand Pivot with Gait Belt	Guarded Stand Pivot Transfer Safety Gait Belt and Baseline Gait Device
			FLOOR	Guarded Stand Pivot Transfer Safety Gait Belt and Baseline Gait Device	Guarded Stand Pivot Transfer Safety Gait Belt and Baseline Gait Device	Guarded Stand Pivot Transfer Safety Gait Belt and Baseline Gait Device
			SUPINE	Resident has demonstrated consistency and independence with the specific transfer task using baseline-related devices.		
			SITTING			
			Standing and Advanced Floor			
	<p>Supine postures include: Bed, Gurney, Diagnostic, and Operative Surfaces</p> <p>Sling Loading Devices include: Bedside Sling Lifts, Wall-mounted Braced Sling Lifts, and A-frame Lifts</p> <p>Slide Glide Devices include: Transfer Sheets, Transfer Mattress, Roller Boards, and Supine Plastic Transfer Boards</p> <p>Sitting Postures include: Wheelchairs, Commodes, Bedside Chair</p> <p>Baseline Devices are those of which the resident has demonstrated consistent use.</p>					